

The BootieButler Shoe Cover System makes dispensing shoe covers a simple, safer and money-saving process

Dispensing shoe covers is now as simple as 1 - 2 - 3:



Place one foot in the dispenser, heel first.



Slowly pull back until cover stretches beyond shoe tip, point toe down and pull back.



Repeat with other foot.

Dispensar las cubiertas de zapatos ahora es tan simple como 1 - 2 - 3:



Coloque un pie en el dispensador, primero el talón.



Lentamente tire hacia atrás hasta que la cubierta se extienda más allá de la punta del zapato, apunte hacia abajo y tire hacia atrás.



Repita con el otro pie.